

Lactate Assay

Biochemistry:

Lactate is oxidized to pyruvate by the lactate dehydrogenase (LDH) reaction. The hydrazine destroys the pyruvate, allowing the reaction to run to the complete oxidation of all lactate molecules. To ensure this, NAD^+ is provided in excess. The concentration of lactate in the sample is proportional to the increase in absorbance as NAD^+ is reduced to NADH.

Sample Preparation:



The blood/tissue sample needs to be deproteinized to prevent the slow release of lactate from erythrocytes. The deproteinization is best done by a small (1:2 or 1:3) dilution in 6% PCA. Alternatively, blood can be treated with an anti-glycolytic agent (sodium fluoride), left to clot, and centrifuged for subsequent collection of serum. The anti-glycolytic agent could also be used with heparin or EDTA to block clotting, with centrifugation and subsequent collection of plasma.

Table 1: Assay cocktail ingredients (designed for glycerol ingestion studies).

Compound	Final []	Amount of Compound for given cocktail volume				
		25 mL	100 mL	150 mL	200 mL	250 mL
Water		16	64	96	128	160
Glycine	320 mmol/L	8 mL	32 mL	48 mL	64 mL	80 mL
Hydrazine	320 mmol/L	0.4 mL	1.6 mL	2.4 mL	3.2 mL	4 mL
NAD^+	2.4 mmol/L	0.00955 g	0.0382 g	0.0573 g	0.0764 g	0.0955 g
Enzymes						
LDH*	2 U/mL	10 μL	40 μL	60 μL	80 μL	100 μL

* LDH stock = 5,435 U/mL

Use distilled water samples for the blanks. Read absorbance at 340 nm, and first zero to the blanks. Calculate lactate from the Beer-Lambert Law, using 6.22 as the millimolar extinction coefficient for NADH. Resting blood samples should approximate 1 mmol/L.

Table 2: Stock Sample Locations

Fridge	Freezer	Bench
Glycine solution (1 M), LDH	NAD^+	Hydrazine (from store room cabinet)